

ARGUS PRO LifeCare 2

Leading-edge technology is the key

Schiller AG surprised experts already in 1974 with an emergency electrocardioscope weighing just 480 grams, which started its triumphal procession throughout the world and laid the cornerstone for countless further successes. Moreover, our 480 grams light pocket defibrillator Easyport is still unrivalled and has saved many lives.

The ARGUS PRO LifeCare 2 heralds a new age in rescue services, meaning that paramedics no longer need to carry unnecessary ballast. Thanks to the skilful use of leading-edge technology and the robust lightweight construction that is also used in aircraft construction, Schiller has designed a defi/monitor with built-in ECG unit that satisfies even highest expectations.

Schiller AG, an autonomous and financially independent company with more than 700 employees and 25 subsidiaries, offers you a worldwide unique service support and guarantees further technical peak performances.

SCHILLER AG pioneering company since 1974

Only 4.4 kg with

8.4" colour display

- Outstanding readability even in sunlight and from a distance
- 6 waveforms and all parameters individually configurable
- Diagnostic ECG with all 12 channels

Soft keys

 7 soft keys for quick and easy operation and direct access to main features


Convenient 6-channel ECG printer

- Paper width: 114 mm
- Paper length: 30 m, Z-folded
- Settings individually configurable
- Paper replacement in a matter of seconds
- Printout is already folded


SPOOL BY THE PARTY OF THE PARTY

Examples of possible display configurations:





all options - makes life easier for paramedics



ARGUS PRO LifeCare 2

12-channel ECG

- Simultaneous 12-channel ECG
- Measurement and interpretation (optional)

NIBP

• Automatic blood pressure measurement for adults, children and neonates

2 USB interfaces

- Device configuration via USB stick
- Software updates
- Memory stick (e.g. for intervention data)
- WLAN

GSM

• Built-in GSM module for data transfer

Connection for defi/pacer pads

• If pads are used, the left paddle connection is covered

Connection for left paddle

• If paddles are used instead of pads



Optimal sto



rage and protection for all accessories



Have everything close at hand in the two side pockets

- Accommodates even the paddles
- All sensors can be connected
- Optimal protection for the plugs



Connections right

ARGUS PRO LifeCare 2

Sophisticated wall mounting

- Extremely light and robust (carbon fibre-reinforced)
- One-hand operation with automatic locking
- Mains connector (100 250 V AC) and 15 Volt DC output
- Connection to on-board power supply 8 48 Volt DC
- Ethernet interface
- RS-232-C interface
- USB interface
- GSM antenna connection
- SIM card slot
- VGA connection for display unit
- Customised mounting on the wall or on horizontal surfaces

If the ARGUS PRO LifeCare 2 is unlocked without being removed from the wall mounting, the unit is again locked after approx. 10 seconds.



Charging unit

The ARGUS PRO LifeCare 2 can be charged without using the wall mounting, by connecting the device directly to the power supply and charging unit, which is supplied as standard. The corresponding connection is situated at the back of the device.







The electrical contacts for the connection with the ARGUS PRO LifeCare 2 are located on the back of the two recording cones and are protected from mechanical damage.



Technical Data

General data:

Dimensions: 260 x 265 x 145 mm (10.2 x 10.4 x 5.7 in)

Weight: 4.4 kg with all options and 2 batteries

Power supply: Lithium ion battery 10.4 V, 4.3 Ah (2x), capacity for at least

100 shocks at 120 J and 6 hours monitoring

Protection class: Applied part monitoring CF, defibrillation and pace-

maker BF

Conformity: MDD 93/42/EEC; class IIb Standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-2-4, IEC 60601-2-27, IEC 60601-2-30, IEC 60601-2-34, IEC 60601-2-49, IEC 9919, IEC 1060-1, IEC 1060-3, IEC 14971, RTCA/DO 160-E, EN 1789, EN 13718

Environmental conditions:

- Temperature, operating: -10° to 40° C, -20° to 60 °C for one hour after device has been stored at room temperature
- Temperature, storage: -10° to 50° C/14° to 122° F
- Relative humidity: 0 to 95 %, non condensing

Display:

LCD colour display 8.4" (211 x 158 mm), transflective, backlit, resolution 1024 x 768 pixels (SVGA), clearly readable even in sunlight

Printing technology: High-resolution thermal printer Chart paper: Thermo-reactive, Z-folded, 114 mm wide

Measured parameters:

Sampling frequency: 1000 Hz

Leads: Simultaneous acquisition of leads I, II, III, aVL, aVF, aVR, V1-V6

Artefact suppression: On/off Pacemaker detection:

- Amplitude: min. 1 to 2 mV
- Atrial and ventricular pacemakers

- Diagnostic ECG:
 Arrhythmia analysis
- Resting ECG including measurement and interpretation
- ST measurement/analysis

Respiration:

Measuring method: Thorax impedance method

- Respiration rate: 0 to 200/min
- Impedance: 0.1 to 3.0 Ohm
- Apnea: 15 to 35 sec

Measuring method: Masimo SET™

Parameters: SpO₂, PP (peripheral pulse) (CO and METHb in preparation)

Parameters: SYS, MEAN, DIA, Pulse

Measurements: Single, automatic, short-time (15 min)

Parameters: SYS, MEAN, DIA, Pulse

Measuring range: SYS, MEAN, DIA: -400 to 400 mmHg

Number of inputs: 2 inputs for simultaneous measurement of 2 tempera-

etCO2 Mainstream and sidestream:

Measuring methods: IRMA™ and ISA™

Range accuracy: CO₂ 0 - 10 % ± 0.3 % ABS or ± 6 %

Respiration rate: ± 1/min

Defibrillator/Transcutaneous pacemaker:

Defibrillation mode: AED and manual mode

Waveform: Multipulse Biowave®

Capacity: 100 shocks at 120J and 6 hours monitoring

AED:

Standard applied energy: 120 - 120 - 120 Joules; different energy levels can be factory-defined

Manual mode:

Energy manual adjustable: 2, 4, 8, 15, 30, 50, 70, 90, 110, 120, 150, 200 Joules

Transcutaneous pacemaker:

Mode: synchronized, non-synchronized Output current: 0 - 150 mA, steps of 5 mA Output rate: 40 - 200/min steps of 5 mA

All specifications are based on an environmental temperature of 25° C

(if not seperately indicated).

Technical data and images are subject to change without notice







SCHILLER Asia-Pacific / Malaysia

52200 Kuala Lumpur, Malaysia Phone +603 6272 3033 Fax +603 6272 2030 sales@schiller.com.my www.schiller-asia.com



Austria

SCHILLER Handelsgesellschaft m.b.H.

A-4040 Linz Phone +43 732 709 90 Fax +43 732 757 000 sales@schiller.at www.schiller.at



SCHILLER Int'l Trading (Shanghai) Co. Ltd.

200335 Shanghai, China Phone +86-21-62099627 Fax + 86-21-62099623 sales@schiller.com.my www.schiller.cn



France SCHILLER Médical S.A.S.

F-67162 Wissembourg/Cedex Phone +33 3 88 63 36 00 Fax +33 3 88 94 12 82 info@schiller.fr www.schiller-medical.com



France (distribution France) SCHILLER France S.A.S.

F-77600 Bussy St Georges Phone +33 1 64 66 50 00 Fax +33 1 64 66 50 10 infoschiller@schiller-france.fr www.schiller-france.com



Germany

SCHILLER Medizintechnik GmbH

D-85521 Ottobrunn Phone +49 89 62 99 81-0 Fax +49 89 609 50 90 info@schillermed.de www.schillermed.de



Hungary SCHILLER Diamed Ltd.

H-1141 Budapest Phone +36 (1) 383-4780 / 460-9491 Fax +36 (1) 383-4778 sales@schiller.at www.schiller-hungary.hu



India

SCHILLER Healthcare India Pvt. Ltd.

Mumbai - 400 001, India Phone +91 22 6152 3333 / 2826 3520 Fax +91 22 2826 3525 sales@schillerindia.com www.schillerindia.com



Japan

SCHILLER Japan, Ltd.

Hiroshima 734-8551 Phone +81 82 250 2055 Fax +81 82 253 1713 koji.maekawa@schiller.jp www.schiller.jp



Poland

SCHILLER Poland Sp. z o.o.

PL-02-729 Warszawa Phone +48 22 8432089 Fax +48 22 8432089 schiller@schiller.pl www.schiller.pl



Russia & C.I.S.

SCHILLER AG Rep. office 125124 Moscow, Russia Phone +7 (495) 970 11 33 Fax +7 (495) 970 11 33 mail@schiller-ag.com

www.schiller-cis.com



Spain SCHILLER ESPAÑA, S.A.

E-28230-Las Rozas/Madrid Phone +34 91 713 01 76 Fax +34 91 355 79 33 schiller@schiller.es www.schiller.es



Switzerland SCHILLER-Reomed AG

CH-8953 Dietikon Phone +41 44 744 30 00 Fax +41 44 740 37 10 sales@schiller-reomed.ch www.schiller-reomed.ch



Turkey SCHILLER TÜRKIYE

Okmeydani-Sisli - Istanbul Phone +90 212 210 8681 (pbx) Fax +90 212 210 8684 sales@schiller-turkiye.com www.schiller-turkiye.com



USA

SCHILLER America Inc.

Doral, Florida 33122 Phone +1 786 845 0620 Fax +1 786 845 06 02 sales@schilleramerica.com www.schilleramerica.com

