


PowerCube[®] ROS



**cooperation-free determination
of the airway resistance**

**as flow/volume test alone is not
sufficient enough!**

GANSHORN
MEDIZIN ELECTRONIC

PowerCube®

The oscillometric resistance-measurement answers the central question for the current degree of obstruction immediately and practically without the patients' cooperation!

A gentle measurement method without influence on the respiratory tract

When measuring with the PowerCube®-ROS »normal« tidal breathing replaces the strenuous, cooperation-dependent, and error-prone breathing manoeuvres. Oscillation volumes smaller than one milliliter make fixing of the cheeks unnecessary.


Compared to the usual methods the bronchial system is practically not being stressed.

Therefore the gentle resistance measurement is ideal for children!



PowerCube®-ROS – how does the measurement work?

Like with the classic method for the determination of airway obstructions the measurement of the oscillatory resistance is based on the physical coherence between pressure (alveolar pressure), flow (breathing flow) and resistance (airway resistance).



For the classic measuring methods the patient produces pressure and flow by his breathing manoeuvre which are being measured and the resistance being calculated from it.

The oscillatory resistance is based on the same physical basis. Here the measuring instrument produces a precisely defined alternating flow (sinus) producing an oscillation pressure in the air ways. The course of oscillation pressure and alternating flow is being measured and the impedance (alternating current resistance) is being calculated from it.

The continuous determination of the phase displacement between oscillation pressure and flow enables diagnostically most interesting conclusions on characteristics of tissue and elasticity of lung and thorax area.





Powerful measuring system for reliable results

- The available measuring frequencies of 5, 10, 20, and 30 Hz are being measured selectively and directly.
The high resolution of up to 60 measurements per second makes the resistance course visible with highest precision as if it was online.
- The detailed resistance volume graphs allow an excellent insight into the breathing mechanics.
- The patient breathes through a reference tube, the defined resistance of which also serves to calibrate the whole system – an important contribution to accuracy of the measurement and to quality assurance.
- The reference impedance guarantees long-time stability and reproducible results since all basis parameters (mouth pressure, flow, and phase displacement) are being measured simultaneously by the automatic calibration.

Determine bronchial obstructions quickly, easily, and surely.

- Cooperation-free test of the bronchial system with reference to resistance, elasticity, and reactivity.
- separation of exobronchial and endobronchial components as alternative to cooperation-depending parameters, e.g. FIV1
- differential diagnostically important information in the form of:
 - resistance-volume graphs,
 - resistance-phase graphs,
 - frequency response of the bronchial system
 - phase response of the bronchial system
- separate resistance determination for inhalation and exhalation to determine bronchial instabilities

Also available as option together with other PowerCube measuring systems.

GANSHORN
MEDIZIN ELECTRONIC

Industriestraße 6-8
D-97618 Niederlauer / Germany
phone: +49 - 9771 - 62 22 - 0
fax: +49 - 9771 - 62 22 - 55
e-mail: sales@ganshorn.de
www.ganshorn.de

PowerCube®-ROS

PC-based measuring module for the test of the bronchial system.

Resistance-measurement including spirometry, flow/volume, maximum voluntary ventilation MVV.

Integrated program for pre-/post-comparisons, bronchospasm tests and provocations

337 020 PowerCube®-ROS with software LF8 and basis accessories
(flow transducer, tubes, mouthpieces, calibration resistors, etc.)

Options

337 105 nose resistance ROS nasal (with face-mask)

223 030 emphysema-diagnosis (CO₂-capnography)

337 150 FRC determination SB and StSt (N₂-wash out)

337 125 rhinomanometry

971 010 movable trolley standard

942 206 med. 230V insulation transformer 1000 VA

The software is network compatible and can be connected by GDT interface with a computer network in the surgery.

870 010 PCS connection

871 040 multi-station license for the establishment of additional examination places

More options, accessories (e.g. trolleys), connections, personal computer requirements on request.

Technical data PowerCube®-ROS:

Active medical product class II a

Dimensions:

Measuring module (depth x width x height): 15 x 15 x 15 cm

Measuring desk (depth x width x height): about. 65 x 120 x 65 cm

Weight:

Measuring module: about 3 kg

Whole measuring desk about 70 kg

Power consumption base unit: about 15 VA

Measuring principle resistance: oscillation (Siemens system)

Measuring ranges: 0.5, 1, 2 und 4 kPa/l/sec.

Oscillation frequencies: 5, 10, 20, 30 Hz

Automatic BTPS correction

RS 232- infra red interface for a 4 KV-insulation between PowerCube and personal computer

Developed, produced and quality-controlled according to
EN ISO 9001/12.2000 / EN ISO 13485/11.2000

CE 0124

Ganshorn Medizin Electronic is not responsible for any errors contained in this document. Images in this document may show optional features at extra charges.

All shown products are subject to technical changes without prior notice. 2/05

