


SCHILLER

DEFIGARD TOUCH 7

Maximum performance, minimum weight,
touch screen control: the power to save lives



SCHILLER
The Art of Diagnostics

DEFIGARD TOUCH 7

A compact combination of state-of-the-art defibrillation technology and comprehensive monitoring functions

The DEFIGARD Touch 7 is an indispensable tool for rescuers. The device is extremely compact and offers the latest defibrillation technology in combination with comprehensive monitoring functions. It is the first emergency monitor/defibrillator equipped with a touch screen – making it the most intuitive device on the market – as well as with the latest data transmission technology.

UNPARALLELED TOUCHSCREEN

The DEFIGARD Touch 7 is the first emergency monitor/defibrillator with a touch screen. Thanks to the new software environment Touch'n Save, this device offers unmatched user-friendliness. Moreover, the exceptional display quality enables optimal display readability at all times. The large 7-inch touch-screen is particularly bright, and its high resolution guarantees unparalleled user-friendliness. Thanks to Projected Capacitive Touch technology, the display enables highly responsive operation, even when wearing gloves.

CAN BE USED ANYWHERE

Thanks to SCHILLER's extensive experience in the rescue sector, the DEFIGARD Touch 7 has been developed so that it can also be used in extreme conditions. It is the ideal companion for all land and water interventions.

Through its lightweight structure, it also offers new possibilities for air rescue operations. Its low weight – three to four times below that of competing products – ensures significant fuel saving when used on an air plane or a helicopter.

With its mounting bracket – which also serves as battery charger – the DEFIGARD Touch 7 is resistant to shocks and vibrations, thus complying with land and air rescue operation standards.

CONNECTED AT ALL TIMES

The device can perform high speed data transmission at all times. Depending on the location, transmission can be performed using Wi-Fi or 3G.

The safety standards WPA2 Enterprise guarantee secure data protection.

PRACTICAL AND SOPHISTICATED

In the functional carrier bag, each accessory has its place, is easily accessible and can readily be used, as each part is already connected. Moreover, the device can be tilted for optimal convenience.

OPERATIONAL AT ALL TIMES

The DEFIGARD Touch 7 is ready to use at all times, thanks to regularly performed self-tests and to a RFID transponder that ensures the monitoring of the defibrillator electrodes' expiry date.





1 Large selection of NIBP cuffs

2 Invasive thermometry or through skin temperature sensor


3 Defibrillation electrodes for children and adults

4 4 + 6 lead ECG cables


5 Masimo Rainbow SpO₂ Technology

COMPACT, LIGHTWEIGHT AND USER-FRIENDLY: YOUR PERFECT COMPANION!

ADAPTIVE USER INTERFACE



PHYSICIAN MODE This display, showing values and waveforms as well as a manual defibrillation function, gives the physician a comprehensive overview of the patient.



RESCUE MODE With a simple swipe gesture, the user can move from the physician display to the rescue display. The values are displayed numerically and are well visible, even from a distance. Moreover, a semi-automatic defibrillation function is also available for all users.

ECG

The DEFIGARD Touch 7 enables continuous ECG monitoring as well as 12-lead ECG diagnostic recording. ECG electrodes are conditioned so as to be used in emergency situations.

RESPIRATION

The respiration rate is measured through impedance variation using the ECG electrodes.

SPO₂, SPCO AND SPMET

The Masimo Rainbow technology for SpO₂, SpCO and SpMet measurements is available as an option.

The carboxyhaemoglobin (SpCO) enables the non-invasive and immediate detection of heightened levels of carbon monoxide in the blood. This ensures a rapid diagnostic and treatment of patients with carbon monoxide poisoning.

The SpMet measurement enables the non-invasive and immediate detection of high levels of methemoglobin in the blood.

NON-INVASIVE BLOOD PRESSURE MEASUREMENT


Arterial blood pressure measurement in children and adults, as well as a venous block function

TEMPERATURE

Continuous measurement of the temperature through skin sensors or endo-cavity probes.




ARGUS PRO LIFEPOINT
The reusable CPR-Feed-back-Sensor is easy to position.



MONITORING
Visual monitoring of the cardiac massage

Manual defibrillation function with cardioversion



ASSISTANCE IN DELIVERING CARDIAC MASSAGE

Thanks to the LifePoint CPR feedback sensor, cardiac massage monitoring is made easy. The following elements are displayed with various colors.

Available information:

- Compression depth and rate
- Pressure relief at the end of each compression¹

MANUAL OR SEMI-AUTOMATIC DEFIBRILLATOR WITH UP TO 200 JOULES

The DEFIBGARD Touch 7 can be operated in semi-automatic or manual mode (synchronously or asynchronously) with Multipulse Biowave (biphasic), which is achieved with exceptional efficacy and safety, as well as with extremely low energy input.

This makes this device ideal for all patient types. It identifies paediatric electrodes and adjusts the emitted energy automatically.

¹ Not valid for Great Britain, Germany and France

DIRECT DATA TRANSMISSION



INTERVENTION SITE

Recording of a 12-channel ECG on site and immediate transmission to the referring physician.




MEDICAL REGULATION

Reading and analysis of the data transmitted from the intervention site.



EMERGENCY TRANSPORT

The emergency doctor can implement the instructions immediately and thus help saving essential minutes.



Direct data transmission ensures optimal care of the patient. The rescuer or physician can remain in contact with the regulating doctor.

The DEFIGARD Touch 7 enables data transmission during and after the intervention.

DURING THE INTERVENTION

- ❖ 12-lead ECG
- ❖ Trends
- ❖ Screenshots
- ❖ PDF reports
- ❖ Connection with ePCR

AFTER THE INTERVENTION

- ❖ Comprehensive intervention data including ECG and SpO₂ curves, events and trends.



CONNECTED AT ALL TIMES

- ❖ This device is compatible with the SCHILLER Update Server, which makes the updating through the network easier than ever before.
- ❖ Manage the entire configuration of the DEFIGARD Touch 7: the devices can be organised into groups, build configuration templates and be deployed easily to many devices at once.
- ❖ Device self-test results are automatically sent by e-mail, for a better pool management and immediate action in case of error.



CONTAINS ALL FUNCTIONS
OF A LARGE MONITORING
DEVICE, BUT WEIGHS
ONLY 2.4 KG



Asia
SCHILLER Asia-Pacific / Malaysia
 52200 Kuala Lumpur, Malaysia
 Phone +603 6272 3033
 Fax +603 6272 2030
sales@schiller.com.my
www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
 A-4040 Linz
 Phone +43 732 709 90
 Fax +43 732 757 000
sales@schiller.at
www.schiller.at



China
Alfred Schiller Medical Equipment Co. Ltd.
 100015 Beijing, China
 Phone +86-10-84565453
info@schillermedical.com
www.schillermedical.cn



France
SCHILLER Médical S.A.S.
 F-67162 Wissembourg/Cedex
 Phone +33 3 88 63 36 00
 Fax +33 3 88 94 12 82
info@schiller.fr
www.schiller-medical.com



France (distribution France)
SCHILLER France S.A.S.
 F-77600 Bussy St Georges
 Phone +33 1 64 66 50 00
 Fax +33 1 64 66 50 10
infoschiller@schiller-france.fr
www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
 D-85622 Feldkirchen b. München
 Phone +49 89 62 99 81 0
 Fax +49 89 609 50 90
info@schillermed.de
www.schillermed.de



Hungary
SCHILLER Diamed Ltd.
 H-1141 Budapest
 Phone +36 (1) 383-4780 / 460-9491
 Fax +36 (1) 383-4778
sales@schiller.at
www.schiller-hungary.hu



India
SCHILLER Healthcare India Pvt. Ltd.
 Mumbai - 400 001, India
 Phone +91 22 6152 3333 / 2920 9141
 Fax +91 22 2920 9142
sales@schillerindia.com
www.schillerindia.com



Croatia
Schiller medicinski instrumenti d.o.o.
 10000 Zagreb
 Phone +385 1 309 66 59
 Fax +385 1 309 66 60
info@schillerz.hr
www.schiller.ch



Latin America
SCHILLER Latin America, Inc.
 Doral, Florida 33172
 Phone +1 954 673 03 58
 Fax +1 786 845 06 02
sales@schillerla.ch
www.schillerla.ch



Poland
SCHILLER Poland Sp. z o.o.
 PL-02-729 Warszawa
 Phone +48 22 843 20 89
 Fax +48 22 843 20 89
schiller@schiller.pl
www.schiller.pl



Russia & C.I.S.
SCHILLER AG Rep. office
 125124 Moscow, Russia
 Phone +7 (495) 970 11 33
 Fax +7 (495) 970 11 33
mail@schiller-ag.com
www.schiller-cis.com



Serbia
SCHILLER d.o.o.
 11010 Beograd
 Phone +381 11 39 79 508
 Fax +381 11 39 79 518
info@schiller.rs
www.schiller.rs



Slovenia
SCHILLER d.o.o.
 2310 Slovenska Bistrica
 Phone +386 2 843 00 56
 Fax +386 2 843 00 57
info@schiller.si
www.schiller.si



Spain
SCHILLER ESPAÑA, S.A.
 E-28230-Las Rozas/Madrid
 Phone +34 91 713 01 76
 Fax +34 91 355 79 33
schiller@schiller.es
www.schiller.es



Switzerland
SCHILLER-Reomed AG
 CH-8953 Dietikon
 Phone +41 44 744 30 00
 Fax +41 44 740 37 10
sales@schiller-reomed.ch
www.schiller-reomed.ch



Turkey
SCHILLER TÜRKİYE
 Okmeydanı-Sisli - İstanbul
 Phone +90 212 210 8681 (pbx)
 Fax +90 212 210 8684
sales@schiller-turkiye.com
www.schiller-turkiye.com



USA
SCHILLER America Inc.
 Doral, Florida 33172
 Phone +1 786 845 06 20
 Fax +1 786 845 06 02
sales@schilleramerica.com
www.schilleramerica.com

