




CARDIOVIT CS-200 Excellence

Electrocardiology of the Future

- 24 bit resolution of the ECG signal
- 8000 Hz sampling frequency
- 16-lead resting and exercise ECG
- Automatic measurement and interpretation
- Intuitive operation
- Unlimited connectivity
- Immediate emergency ECG



SCHILLER

The Art of Diagnostics

Electrocardiology of the Future

SCHILLER CARDIOVIT® CS-200 Excellence

SCHILLER, in co-operation with world-leading cardiologists and physicians, has developed for you an innovative leading-edge ECG diagnostic system. This is the latest innovation of SCHILLER, the pioneer in computerised ECG measurement and interpretation with more than 35 years of experience in the ECG domain.

The CARDIOVIT® CS-200 Excellence sets new standards for the future in electrocardiology!

- Intuitive and user-friendly operation
- 24-bit amplitude resolution with 8000 Hz sampling frequency
- 12/16 channels for resting and exercise ECG

- Autonomous emergency ECG (unique world-wide)
- SCHILLER ETM® analysis – automatic ECG measurement and interpretation
- Online rhythm ECG for 12/16 channels
- SCHILLER ALPS® – automatic alignment and coloured visual representation of all heart beats during exercise
- Full integration with practice and clinic systems
- Flexible data export and import




Leading-Edge Technology

SCHILLER sets new standards in technology and know-how.



Connectivity

Ethernet and wireless functions enable the bi-directional communication with practice and clinic systems (e.g. via the SEMA® cardiology information system). Patient data can very quickly and easily be viewed, managed and archived.



Top Quality

Perfect presentation and evaluation of your resting and exercise ECG signals in up to 16 channels.

Innovation for highest demands!

Resting ECG

The CARDIOVIT® CS-200 Excellence is the ideal system for extraordinary 12/16-lead resting ECG recordings.

Its intuitive user interface is the basis for highest productivity and efficiency.

High Signal Quality

An extremely high amplitude resolution combined with a high sampling frequency ensure perfect ECG recordings.

Ready for Use at any Time

The integrated thermal printer and rechargeable battery guarantee immediate real-time printouts at any time (autonomous emergency ECG).

Open for the Data World

The CARDIOVIT® CS-200 Excellence can be integrated with any network. Central archiving and data access is possible from anywhere.

Automatic Measurement and Interpretation

The proven SCHILLER measurement and interpretation analysis supports you with interpretation statements for the evaluation of resting ECGs – for adults and children.




Highest Signal Quality

The automatic signal quality monitoring function already informs you about any disturbances while you are applying the electrodes. You have full control over the signal quality.



Printout on Thermal or Normal Paper

You select if you want to print 3, 6, 12 or 16 channels on thermal or normal paper – in top quality.



New ECG Analysis Standard

The SCHILLER ETM® ECG analysis program provides you with all the information you need.



Exemplary Ergometry Solutions

Exercise ECG

The CARDIOVIT® CS-200 Excellence ensures the continuous monitoring of your patient during the entire exercise ECG.

Safe Diagnosis

The extremely high amplitude resolution and sampling frequency allow the very accurate detection of ST alterations and Duke treadmill score calculation. Precise blood-pressure values ensure the correct diagnosis.

Optimised Work Process

You can already assess an ECG and write the interpretation during the exercise.

Online Exercise ECG

On your demand, all 12/16 leads over the entire ergometry will be saved.


SCHILLER ALPS®

The SCHILLER ALPS® gives you the visual presentation of QT intervals, PR intervals and ST alteration at any time during an exercise ECG



Function Keys

The direct function keys allow a swift ergometry procedure. They provide the convenient and straightforward control of the exercise ECG.



Continuous Monitoring

Continuous average updates, ST amplitudes, ST/HR slope, ST/HR loops and ST/HR index support your diagnosis.



ST alterations

With the SCHILLER ALPS®, you can immediately detect ST alterations at any time during an exercise ECG.



Asia

SCHILLER Asia-Pacific / Malaysia
52200 Kuala Lumpur, Malaysia
Phone +603 6272 3033
Fax +603 6272 2030
sales@schiller.com.my
www.schiller-asia.com



Austria

SCHILLER Handelsgesellschaft m.b.H.
A-4040 Linz
Phone +43 732 709 90
Fax +43 732 757 000
sales@schiller.at
www.schiller.at



China

SCHILLER Int'l Trading (Shanghai) Co. Ltd.
200335 Shanghai, China
Phone +86-21-62099627
Fax + 86-21-62099623
sales@schiller.com.my
www.schiller.cn



France

SCHILLER Médical S.A.S.
F-67162 Wissembourg/Cedex
Phone +33 3 88 63 36 00
Fax +33 3 88 94 12 82
info@schiller.fr
www.schiller-medical.com



France (distribution France)

SCHILLER France S.A.S.
F-77600 Bussy St Georges
Phone +33 1 64 66 50 00
Fax +33 1 64 66 50 10
infoschiller@schiller-france.fr
www.schiller-france.com



Germany

SCHILLER Medizintechnik GmbH
D-85521 Ottobrunn
Phone +49 89 62 99 81-0
Fax +49 89 609 50 90
info@schillermed.de
www.schillermed.de



Hungary

SCHILLER Diamed Ltd.
H-1141 Budapest
Phone +36 (1) 383-4780 / 460-9491
Fax +36 (1) 383-4778
sales@schiller.at
www.schiller-hungary.hu



India

SCHILLER Healthcare India Pvt. Ltd.
Mumbai - 400 001, India
Phone +91 22 6692 0520/ 2826 3520
Fax +91 22 2826 3525
sales@schillerindia.com
www.schillerindia.com



Japan

SCHILLER Japan, Ltd.
Hiroshima 734-8551
Phone +81 82 250 2055
Fax +81 82 253 1713
koji.maekawa@schiller.jp
www.schiller.jp



Poland

SCHILLER Poland Sp. z o.o.
PL-02-729 Warszawa
Phone +48 22 8432089
Fax +48 22 8432089
schiller@schiller.pl
www.schiller.pl



Russia & C.I.S.

SCHILLER AG Rep. office
125124 Moscow, Russia
Phone +7 (495) 970 11 33
Fax +7 (495) 970 11 33
mail@schiller-ag.com
www.schiller-cis.com



Spain

SCHILLER ESPAÑA, S.A.
E-28230-Las Rozas/Madrid
Phone +34 91 713 01 76
Fax +34 91 355 79 33
schiller@schiller.es
www.schiller.es



Switzerland

SCHILLER-Reomed AG
CH-8953 Dietikon
Phone +41 44 744 30 00
Fax +41 44 740 37 10
sales@schiller-reomed.ch
www.schiller-reomed.ch




Turkey

SCHILLER TÜRKİYE
Okmeydanı-Sisli - İstanbul
Phone +90 212 210 8681 (pbx)
Fax +90 212 210 8684
sales@schiller-turkiye.com
www.schiller-turkiye.com



USA

SCHILLER America Inc.
Doral, Florida 33122
Phone +1 786 845 0620
Fax +1 786 845 06 02
sales@schilleramerica.com
www.schilleramerica.com



SCHILLER

The Art of Diagnostics

