

Phyaction E

Phyaction E, for electrotherapy.

With the Phyaction series you treat your patients in a simple and effective way.

- Large and crystal-clear LCD screen
- Knob for precise intensity adjustment
- Simple and easy to use
- Extensive indication and objectives menu
- Compact unit and unique design



Article number: 323 254

Characteristics

- Therapies :
 - o Electro therapy, 2 and 4 poles, 2 channels completely independent
 - o Simultaneous therapy (2 different indications treated simultaneously by using electro)
- Currents :
 - o Current forms: see on the next page
- Functionalities :
 - o Objectives: 63
 - o Indication: 134
 - o Diagnostic: 9
 - o Free memory: 50
- Phyaction E can be connected to Phyaction V

Technical specifications

Dimensions Phyaction	: 265 x 275 x 122 mm
Weight Phyaction	: 3.650 kg
Weight including accessories	: 4.6 kg
Mains voltage	: 100 - 240 VAC, 50-60 Hz
Maximum power, in operation	: 85 VA
Safety class	: Class I (earthed socket required)
Insulation	: Type BF (floating patient circuit)
Fuses	: 2 x T2AL250 V
Treatment time	: 0 - 60 min.
Languages	: English, Dutch, German, French, Italian, Spanish and Portuguese

Standard accessories

108 725	Electrode cables
109 959	Electrode rubber 6 x 8 cm
100 658	Chamex bag for electrode 6 x 8 cm (per 4)
108 935	Fixing strap, elastic 5 x 60 cm
100 689	Mains lead
108 919	Test plug F/F 4 mm
115 684	VAS

Manual

322 890	Dutch
322 934	French
322 846	English
322 978	German
327 415	Spanish
323 000	Safety instruction

Optional accessories

325 068	BAG (carrying case)
329 978	Vaginal probe Novatys Gold
330 594	Vaginal probe V2B+
330 572	Vaginal probe Optima 3
330 583	Vaginal probe Perisize 4+
329 989	Anal probe Analia
330 561	Anal probe Analys+
112 166	Rectal probe
109 958	Electrode rubber 4 x 6 cm (per 2)
109 960	Electrode rubber 8 x 12 cm (per 2)
100 657	Chamex bag for electrode 4 x 6 cm (per 4)
100 659	Chamex bag for electrode 8 x 12 cm (per 4)
326 799	Electrode, adhesive, reusable - ø 3 cm (per 4)
326 810	Electrode, adhesive, reusable - 2,5 x 5 cm (per 4)
326 821	Electrode, adhesive, reusable - 5 x 5 cm (per 4)
326 832	Electrode, adhesive, reusable - 5 x 10 cm (per 4)
328 504	Adhesive electrode kit
113 334	Adapter for adhesive electrode 4 > 2 mm
114 142	Pen electrode with chamex bag
109 944	Chamex for pen electrode 10 pieces
108 934	Fixing strap 5 x 30 cm
108 936	Fixing strap 5 x 120 cm
117 210	Mohil 700





Phyaction E : Current forms
















Unidirectional Currents		
Medium frequency rectangular	Rectangular pulse	Triangular pulse
2-5 current (Ultra Reiz)		




Iontophoresis	
Iontophoresis medium frequency rectangular current	

Diagnostic programs
S/d curve

Diadynamic Currents		
		
MF	DF	CP
		
CP-ISO	LP	LP-ISO

Tens Currents		
		
Conventional Tens example: Symmetrical	Low Frequency Tens example: Symmetrical	Random frequency Tens example: symmetrical
		
Burst Tens example: symmetrical	Han Stim	

NMES Currents		
		
Rectangular Surge	Triangular Surge	Biphasic Surge example: symmetrical
		
Biphasic Surge example: a-symmetrical	Intrapulse Interval Surge Symmetrical	Russian Stimulation
		
2-pole Medium Frequency Surge	Classical Interferential surge	Isoplanar Vector field surge

Interferential Currents		
		
2-pole Medium frequency	Dipole vectorfield	Isoplanar Vectorfield