
Height-adjustable verticalisation table with two motors. One motor ensures effortless and smooth tilting; the other adjusts height and facilitates patient transfer. The Tilting is the ideal standing aid for neuromotor revalidation and adaptation exercises from lying to standing. Another area where it can be applied is general orthostatic examination, more specifically examination and treatment of vasovagal syncope.

- Motors are operated by a hand switch.
- Two separately adjustable fold-away footrests.
- Individually lockable large casters guarantee the best possible mobile use.
- Fitted with a tilt angle indicator.
- Suitable for use with a patient hoist.
- Optional: battery power; separate leg fixation.











## Standard features

- Cushions in a choice of 8 colours.
- Cushions with padding of 40 mm.
- Suitable for use with a patient hoist (clearance below frame = 14cm).
- Gymna Stability Profile (GSP).



## Ordering data - Tilting Advanced



121 095 Electrical with handswitch

Specifications

Electric motor strength : 8.000 N Power supply : 110 - 240 VAC

: 110 - 240 VAC 50/60 Hz or 100 VAC 50/60 Hz

Maximum consumption : 0,8 A Insulation class : II Moister resistance : IP 54 Classification : I

Conformity : CE according to the Medical

Devices Directive of the European Commission

93/42/EEG

Guarantee : 3 Years (upholstery 1 year)

## Optional

125 606 Separate leg fixation125 607 back-up battery power



## Accessories

112 218 Work shelf for tilting125 655 Armrests (set of 2)112 215 Leg fastening strap, large

115 768 Leg fastening strap, for separate leg sections (set of 2)

112 214 Pelvic/thoracic fastening strap

Terry-cloth covers (9 colours available)