

### **Traction** Classic


Special treatment couch for lumbar and cervical traction therapy with a traction device. The large sliding cushion eliminates friction between the patient and the couch surface, resulting in an efficient and comfortable treatment. The support surface can also easily be fixed in position, so that the couch can be used for other purposes as well.




- Head support can be adjusted positively and negatively
- Suitable for use with a patient hoist. (clearance below frame = 15cm).



- Standard features
- All sections gas spring supported.
- Wheel-lifting mechanism with 2 casters and fixed wheels in the front.
- Cushions in a choice of 8 colours.
- Cushions with padding of 40 mm.
- Suitable for patient hoists up to 15cm.



	A La	
A L		J
and the second		
(C)		




lersion number: 1.2

## Autotrac 460 Traction unit

The Autotrac 460 is a digital traction unit for continuous, intermittent or harmonic intermittent lumbar and cervical traction-therapy. It is a nononsense machine with all must-have functionality and rock solid build quality. Together with the great user friendliness, high safety standards and nice design the Autotrac 460 will meet your highest professional requirements.





#### Standard features

- Three therapy modes: continuous, intermittent and harmonic intermittent traction.
- Automatic tension compensation.
- Maximum traction force adjustable in three levels: 1 - 50 or 91 kg for lumbar traction, and 1 - 18 kg for cervical traction.
- Stepless adjustable traction speed.
- Timer from 0 -99 minutes, with remaining time displayed.
- Traction force displayed in kg, N or Ib.
- Digital display for preset force and actual force.
- Digital displays for treatment phase (progressive traction - regressive - pause).
- Buzzer with adjustable sound level for end of treatment or protection alarms.

#### Ordering data

 333.245 Autotrac 460 Traction unit Standard accessories: Mains lead (straight, Europe) Emergency stop switch Extension cord with hooks Four knobs for fixation to the Extension frame Treatment Guide User manual

## 

specifications	
Dimensions (w x I x h) :	
Weight :	
Mains voltage :	
Power consumption :	
Insulation class :	
EMC tested :	
Classification: :	
CE marking :	

# ± 30 x 36 x27 cm ± 13.6 kg 220 - 240 VAC 50/60 Hz Max. 150 VA I type BF (IEC 60601-1) IEC 60601-1-2 Class IIb MDD 93/42/EEC Annex IX CE 0197

#### Optional accessories

100 457	neck harness according to Glisson
100 458	Transverse suspension bar
100 459	Set of washable pads
100 476	Hip belt with Velcro fastening
100 477	Chest belt with Velcro fastening
110 487	Flexion stool
121 434	Extension frame for Autotrac 450/460



